




# Technology

The technique involves the injection of certain substances, via a very fine needle, underneath the dermal layer of the skin into the subcutaneous fat tissue in the area selected. The combination of substances injected depends on the concoction prepared by the practitioner or 'mesotherapist'. It may contain minute quantities of plant extracts, vitamins, enzymes, nutrients, hormones and medications such as vasodilators, non-steroidal anti-inflammatory agents, beta-receptor agonists and antibiotics. Claims of the benefits of dermal injections of hyaluronic acid, antioxidants, dimethylaminoethanol, procollagen peptide, silicon, minerals, coenzymes, amino acids and multi-vitamins include improved elasticity, collagen synthesis, radiance, circulation, homeostasis, micro-exfoliation, skin thickness, hydration and skin softness. The technique involves the injection of certain substances, via a very fine needle, underneath the dermal layer of the skin into the subcutaneous fat tissue in the area selected. The combination of substances injected depends on the concoction prepared by the practitioner or 'mesotherapist'. It may contain minute quantities of plant extracts, vitamins, enzymes, nutrients, hormones and medications such as vasodilators, non-steroidal anti-inflammatory agents, beta-receptor agonists and antibiotics. Claims of the benefits of dermal injections of hyaluronic acid, antioxidants, dimethylaminoethanol, procollagen peptide, silicon, minerals, coenzymes, amino acids and multi-vitamins include improved elasticity, collagen synthesis, radiance, circulation, homeostasis, micro-exfoliation, skin thickness, hydration and skin softness.




### Handpiece



#### surface, so that the skin pores can be quickly opened, and the skin beauty products quickly penetrate into the skin cell layer following the channel. NCIPLE intained

## Functionality





## Benefit



RS rejuvenation-01		
Parameters	Quantity	
Туре	Water Mesotherapy Gun	
Multi-needle	5needles, 9needles	
Mode	Low speed/ High speed/ 10-250 times/min	
Injector capacity	1ml, 2ml, 5ml	
Advantage	Single needle, multi needle, vacuum, 3 in 1	
Vacuum flow	25L/min pressure:90Kpa	
Two modes	Single mode and continuous mode	
Depth	0.3-4.5mm	
Vacuum	1-4pa	
Speed	1-6s	
Dose	4-64μΙ	
Display	4.3 inch LCD touch screen	
Power	50W	
Voltage	220V, 50Hz	
DA	12V/5A	